

## Avian DNA Tests Submission Form

Contact Information	
<b>Owner/Breeder</b>	<b>Veterinarian</b>
Name:	Dr.:
Company:	Clinic:
Address:	Address:
City:	City:
State/Province:	State/Province:
ZIP/Postal Code:	ZIP/Postal Code:
Country:	Country:
<b>Report by:</b>	
<input type="checkbox"/> Phone: _____	<input type="checkbox"/> Fax: _____
<input type="checkbox"/> E-mail: _____	<input type="checkbox"/> Other: _____

Avian DNA Tests					
Code	Test	Sample	Code	Test	Sample
<b>A101</b>	Avian DNA Sexing	Blood, Cards, Feathers	<b>A118</b>	<i>Avian Tuberculosis</i>	Blood, Tissue, Feces
<b>A102</b>	Avian DNA Fingerprinting / Family Relation Test	Blood, Feathers	<b>D409</b>	Aspergillus spp.	Blood, Tissue
<b>A104</b>	Avian Polyomavirus	Blood, Tissue, Feces	<b>D411</b>	West Nile virus	Blood, Tissue, Swab
<b>A1064</b>	<i>Chlamydophila psittaci</i> (Psittacosis)	Blood, Tissue, Feces, Swab	<b>D417</b>	<i>Bordetella avium</i>	Blood, Swab
<b>A108</b>	Psittacine Beak and Feather Disease (Pbfd) virus	Blood, Tissue, Feces	<b>FAMILY PROTECTION Health Program</b>		
<b>A110</b>	Pacheco's Disease virus	Blood, Tissue, Feces	<b>AZP2</b>	Avian Zoonotic Profile: <i>Chlamydophilal psittaci</i> , <i>Cryptosporidium</i> spp., Avian Tuberculosis ( <i>Mycobacterium</i> spp.), <i>Toxoplasma gondii</i>	Blood, Feces
<b>A114</b>	<i>Trichomonas gallinae</i>	Feces			
<b>A116</b>	Avian Paramyxovirus	Blood, Tissue, Feces			

Sample Information	
<p>1. Bird ID: _____ Species: _____ Owner: _____</p> <p><input type="radio"/> A101   <input type="radio"/> A102   <input type="radio"/> A104   <input type="radio"/> A1064   <input type="radio"/> A108   <input type="radio"/> A110   <input type="radio"/> A114 <input type="radio"/> A116   <input type="radio"/> A118   <input type="radio"/> D409   <input type="radio"/> D411   <input type="radio"/> D417   <input type="radio"/> AZP2   <input type="radio"/> ____</p>	<p>2. Bird ID: _____ Species: _____ Owner: _____</p> <p><input type="radio"/> A101   <input type="radio"/> A102   <input type="radio"/> A104   <input type="radio"/> A1064   <input type="radio"/> A108   <input type="radio"/> A110   <input type="radio"/> A114 <input type="radio"/> A116   <input type="radio"/> A118   <input type="radio"/> D409   <input type="radio"/> D411   <input type="radio"/> D417   <input type="radio"/> AZP2   <input type="radio"/> ____</p>
<p>3. Bird ID: _____ Species: _____ Owner: _____</p> <p><input type="radio"/> A101   <input type="radio"/> A102   <input type="radio"/> A104   <input type="radio"/> A1064   <input type="radio"/> A108   <input type="radio"/> A110   <input type="radio"/> A114 <input type="radio"/> A116   <input type="radio"/> A118   <input type="radio"/> D409   <input type="radio"/> D411   <input type="radio"/> D417   <input type="radio"/> AZP2   <input type="radio"/> ____</p>	<p>4. Bird ID: _____ Species: _____ Owner: _____</p> <p><input type="radio"/> A101   <input type="radio"/> A102   <input type="radio"/> A104   <input type="radio"/> A1064   <input type="radio"/> A108   <input type="radio"/> A110   <input type="radio"/> A114 <input type="radio"/> A116   <input type="radio"/> A118   <input type="radio"/> D409   <input type="radio"/> D411   <input type="radio"/> D417   <input type="radio"/> AZP2   <input type="radio"/> ____</p>
<p>5. Bird ID: _____ Species: _____ Owner: _____</p> <p><input type="radio"/> A101   <input type="radio"/> A102   <input type="radio"/> A104   <input type="radio"/> A1064   <input type="radio"/> A108   <input type="radio"/> A110   <input type="radio"/> A114 <input type="radio"/> A116   <input type="radio"/> A118   <input type="radio"/> D409   <input type="radio"/> D411   <input type="radio"/> D417   <input type="radio"/> AZP2   <input type="radio"/> ____</p>	<p>6. Bird ID: _____ Species: _____ Owner: _____</p> <p><input type="radio"/> A101   <input type="radio"/> A102   <input type="radio"/> A104   <input type="radio"/> A1064   <input type="radio"/> A108   <input type="radio"/> A110   <input type="radio"/> A114 <input type="radio"/> A116   <input type="radio"/> A118   <input type="radio"/> D409   <input type="radio"/> D411   <input type="radio"/> D417   <input type="radio"/> AZP2   <input type="radio"/> ____</p>
<p>7. Bird ID: _____ Species: _____ Owner: _____</p> <p><input type="radio"/> A101   <input type="radio"/> A102   <input type="radio"/> A104   <input type="radio"/> A1064   <input type="radio"/> A108   <input type="radio"/> A110   <input type="radio"/> A114 <input type="radio"/> A116   <input type="radio"/> A118   <input type="radio"/> D409   <input type="radio"/> D411   <input type="radio"/> D417   <input type="radio"/> AZP2   <input type="radio"/> ____</p>	<p>8. Bird ID: _____ Species: _____ Owner: _____</p> <p><input type="radio"/> A101   <input type="radio"/> A102   <input type="radio"/> A104   <input type="radio"/> A1064   <input type="radio"/> A108   <input type="radio"/> A110   <input type="radio"/> A114 <input type="radio"/> A116   <input type="radio"/> A118   <input type="radio"/> D409   <input type="radio"/> D411   <input type="radio"/> D417   <input type="radio"/> AZP2   <input type="radio"/> ____</p>

<p><b>Notes / Comments / Requests:</b></p>	<p><b>Payment:</b></p> <p>Credit Card #: _____</p> <p>Exp. Date: _____ Amount: _____</p> <p>Signature: _____ HST (13%) - Canada only: _____</p> <p>Total Amount: _____</p> <p style="text-align: right;"> <input type="checkbox"/> MasterCard  <input type="checkbox"/> Cheque  <input type="checkbox"/> Cash  <input type="checkbox"/> VISA  <input type="checkbox"/> AmEx         </p>
--	--

## **Blood Collection Instructions:**

1. Write the Bird's ID and your name on the card provided. Since the toenail is the easiest location to collect blood, wipe the bird's nail with an alcohol swab and let dry (5 - 10 seconds).
2. Clip the nail short enough to allow bleeding (avoid clipping too high since this may be uncomfortable for the bird). If blood is not present, gently squeeze the toe until blood appears. Ensure that the toenail clipper is thoroughly cleaned with alcohol between birds.
3. Allow a small amount of blood to bead up and carefully touch the corner of the card to the blood bead (only one or two drops is required). Sample only one bird per card.
4. Leave the card to dry for approximately 20 minutes at room temperature. Once the card is dry, place it into the small envelope provided (only 1 card per envelope).
5. Once the samples have dried and the cards have been placed in their individual envelopes and sealed, the envelopes may be collected and placed into a mailing envelope. Ensure ID information is on the cards and the submission form provided. If you are not submitting the sample(s) immediately, store the cards at room temperature in a dark dry location. Cards can be stored in this fashion for several weeks.

**Price: \$20.00 CDN + 13% HST**

## **Feather Collection Instructions:**

1. Select only medium sized mature feathers for plucking. The feathers must be physically plucked (molted feathers cannot be used for this test). Very small feathers located under the wing are not recommended.
2. Pluck 4-6 medium sized feathers from the chest area of the bird and place them in a zip-lock bag. The quantity of feathers plucked is not as important as the quality. Ensure that the feather tip, which is attached to the skin, is present (white colouration). Large feathers may contain some blood which should be left on the feather.
3. Sample only 1 bird's feathers per zip-lock bag. Mark the bird's identification number or name directly on the bag with magic marker (alternatively, place a piece of paper with the birds name and species in the bag). Ensure that the submission form is properly filled out including your name, address (and other contact information), as well as the bird's name and species.

**Price: \$24.00 CDN + 13% HST**