

Heartworm Ag Test Kit Order Form

Canine Heartworm Ag Rapid Test Kit is a sandwich lateral flow immunochromatographic assay for the qualitative detection of *Dirofilaria immitis* in dog's blood samples. Please review the Heartworm Ag Rapid Test Kit specification below:

1. Assay time: **5 min.**
2. Specificity: 99.8%.
3. Sensitivity: 93.5%.
4. Sample: serum, plasma or whole blood.
5. Shelf life: **18 months under room temperature.**
6. ISO 9001, 13485 certified.

Simply complete this Heartworm Ag Test Kit Order form and fax to **289-553-5232**.

Clinic: _____ Account Number: _____

Contact Person: _____ Title: _____

Address: _____

City: _____ Province: _____ Postal Code: _____

Phone: _____ Fax: _____

Please choose one of the following options:
(minimum order is 1 box [**10 kits per box**])

Term	Total # of Boxes
<input type="checkbox"/> Orders of 1 - 4 boxes	_____
<input type="checkbox"/> Orders of 5 - 19 boxes	_____
<input type="checkbox"/> Orders of 20 boxes or more	_____

Payment in Full: Add to the clinic invoice

Cheque enclosed

Credit Card:

_____ Exp. _____

Signature: _____ Date: _____